



ClearScan

Radio-translucent electrodes for ECG monitoring in CathLab, X-Ray and thallium testing procedures


Keep your imaging free of obstruction by using Vermed's ClearScan electrodes. Designed with a carbon-filled plastic snap that will not interfere with your Nuclear Imaging Systems.

ClearScan electrodes are available in foam, cloth or tape, and in a variety of sizes and packaging configurations.

Vermed also carries a complete line of radio-translucent leadwires, as well as a complete line of recording chart paper.

Ask your representative for more information.

- **Radio translucent Ag/AgCl sensor**
- **Hypoallergenic adhesive**
- **Latex-free**
- **Available in economically packaged, no-waste quantities**



Part Number	Description	Number of Electrodes			
		Per Strip	Per Pouch	Per Box	Per Case
Wet Gel					
A10004-5	1 1/2" Diameter Foam	5	5	100	600
A10004-60	1 1/2" Diameter Foam	5	60	-	600
A10004-5T	1 5/8" Diameter Tape	5	5	100	600
A10004-60T	1 5/8" Diameter Tape	5	60	-	600
A10006-5RT	1 7/8" Diameter Foam	5	5	100	600
A10006-1-60RT	1 7/8" Diameter Foam	1	60	-	600
A10025-60	1 7/8" Diameter Foam	5	60	-	600
A10025T-60	1 7/8" Diameter Tape	5	60	-	600
Solid Gel					
A10031-60RT	1 7/8" Diameter Foam	5	60	-	600
A10032-60RT	1 7/8" Diameter Tape	1	60	-	600

To place an order or receive sample electrodes call your sales representative or Customer Service today!

Tenby Medical

A DIVISION OF PEKISKO HOLDINGS INC.

443 Village Square
Sherwood Park, AB
Canada T8A 4K2

Tel/Fax: 780-570-5755
tenby-medical.com