



**Media Archival and Management Systems**

**Turnkey Solution • Store over 250 Hours of Video and Thousands of Images**

**Works With All Medical Cameras • State of the Art Image Display • Side-By-Side Comparison**

**One-Step Media Report • Easy Editing for Presentations • FDA Class 1 Certified**

**Our new D-Scope® System Version 5 family of products includes:**

Capture Workstation, Archive Server and Server Software, Web Browser and iPad Application, HD Acquisition and Streaming.

**With new and updated features including:**

Customized Reporting, Advanced Search, Graphical and Data Annotation, Presentation Editing.

**The D-Scope® System Product Family**

The D-Scope® System Version 5 is a scalable family of Electronic Image Archive products that provide everything needed to capture and archive video and images, and much more.

**Features**

**Capture**

- Full screen video and image capture
- Foot pedal, iPad and voice control
- Standard and high definition video
- Image and Video format flexibility

**Store**

- Locally on a stand-alone system
- Centrally on a variety of storage solution

**Outcome Tools**

- View multiple media files side by side
- Attach diagnostic keywords and comments
- Customizable reporting templates
- Advance search across patients
- Graphical and data annotation tools

**Import and Export**

- Import: DICOM Radiology, image, and digital video files from operating room or other sources
- Presentation Editing with Export to: iPad, Web, video editor, PowerPoint®, CD/DVD.

**D-Scope Workstation**



Our stand-alone D-Scope® Workstation is ideal for private practice, institutional and Operating Room use.

The D-Scope® Workstation is a reliable, high performance, turn-key video/image recording and archival system. The system includes; a powerful medical grade PC with fault tolerant media storage configurable with up to 2 TB of on-line storage; a DICOM compliant 19" or optional 24" monitor; and a durable video cart designed to stand up to high usage in a medical environment.

Our storage solutions provide the ability to configure multiple Workstations on a local network, and to share a common patient database on a central media archive.



### Service and Support

We pride ourselves on being client-service driven. Our systems come with 1 year of support and software upgrades including full parts replacement. We offer 3 year extended support renewals and affordable hardware upgrades after 3 years to extend the life of your purchase.

### Scalable Storage Solutions:

D-Scope® **Media Management System (MMS)** A PC Desktop Server with D-Scope® MMS Software. RAID storage of up to 2 Terabyte.

D-Scope® **Enterprise Storage System (ESS)** A scalable storage solution with enterprise class hardware redundancy and 7 x 24 hour support with D-Scope® ESS Software. Storage capability begins at 5 Terabyte.

D-Scope® **Enterprise Storage System (ESS) Software-Only** Our ESS Software installed on your institution's enterprise-class server or VMware.

## Integrate the D-Scope® with any EMR System

### EMR/EHR Integration Modules

D-Scope® **ADT Patient and Scheduling Module** Synchronize patient account and scheduling information between the D-Scope® and an EMR system or ADT Stream. Display patient look up from EMR with Schedule or patient views including MRN.

D-Scope® **2-Way Transfer Module** Integrate media and diagnostic information acquired with a D-Scope® Workstation into the Electronic Medical Record.

Transfer patient demographic information from the EMR to the D-Scope® for inclusion in diagnostic/procedure reports.

**Caution:** Federal law restricts this device to sale by or on the order of a physician.

**MDDDev**  
www.mddev.com

**Medical Digital Developers, LLC**  
(toll free) 877-223-2660

65 9th St, Suite 13-2,  
Brooklyn, NY 11215  
718-768-2660

Distributed by:

**Tenby Medical**

A DIVISION OF PEKISKO HOLDINGS INC.

443 Village Square  
Sherwood Park, AB  
Canada T8A 4K2

Tel/Fax: 780-570-5755  
tenby-medical.com