



Vermed EMS Electrodes


When every second counts...

Vermed EMS electrodes are designed and packaged specifically for Emergency Services. Large size and aggressive adhesive insure quick placement and fast, accurate readings. All Vermed electrodes use the highest quality materials and meet or exceed AAMI specifications for disposable ECG electrodes.

- **2 1/4" diameter provides large adhesive area**
- **High chloride wet gel gives fast, clear traces**
- **Can be defibrillated**
- **Packaged in convenient, no-waste quantities for EMS use**
- **Hypoallergenic and Latex Free**
- **Electrodes available in other sizes and materials**



New Packaging Provides Instant Identification!



Specially designed EMS 10-Pack



Order Number	Description	Per Strip	Number of Electrodes		
			Per Pouch	Per Box	Per Case
A10003-3	2 1/4" Foam, Wet Gel	3	3	75	600
A10003-60	2 1/4" Foam, Wet Gel	3	60	–	600
A10003-4	2 1/4" Foam, Wet Gel	4	4	100	600
A10003-4-60	2 1/4" Foam, Wet Gel	4	60	–	600
A10005-4NB	1 1/2" Foam, Wet Gel	4	4	–	600
A10012-3SNB	Large Rectangular Foam, Wet Gel	3	3	–	600
A10012-6S	Large Rectangular Foam, Wet Gel	6	6	–	600
A10052-10S	EMS Ten Pack for 12 Lead ECG	50 Kits per Case, 500 electrodes			

To place an order or for more information, call your sales representative or customer service today!

Tenby Medical

A DIVISION OF PEKISKO HOLDINGS INC.

443 Village Square
Sherwood Park, AB
Canada T8A 4K2

Tel/Fax: 780-570-5755

tenby-medical.com