

ENDODIGI DIGITAL DATA DOCUMENTATION SYSTEM

PICTURE / AUDIO / VIDEO



- PORTABLE
- USER FRIENDLY

ENDODIGI | DIGITAL DATA DOCUMENTATION SYSTEM



- **PORTABLE**
- **USER FRIENDLY**
- **EXTERNAL USB**
- **FOOTSWITCH**
- **INPUT / OUTPUT**
- **PATIENT DATA BASE**
- **PATIENT REPORTS WITH PICTURES**
- **AUDIO AND VIDEO FILES**



MINIMUM REQUIREMENTS FOR RUNNING ENDODIGI SOFTWARE

- OPERATING SYSTEM: - WINDOWS XP OR WINDOWS VISTA
- PROCESSOR: - PENTIUM IV 2.2 GHZ OR FASTER
- ATHLON 64 2.0 GHZ OR FASTER
- HARD DRIVE: - 80Ghz (or LARGER)
- PORTS USB 2.0: - 2 AVAILABLE
- SCREEN RESOLUTION: - 1024 X 768
- CD-RW: - RECOMMENDED FOR BACKUP

TECHNICAL SPECIFICATIONS

MODEL	ENDODIGI
CLASS	1
WEIGHT	1.5 Kg 3 Lb.
DIMENSIONS	12 x 4 x 25 cm. 4.72 x 1.57 x 9.84"
NORMS	PAL NTSC
RESOLUTION	24 BITS (RGB)
FILE FORMATS	JPEG (PICTURES) WMV9 MPEG4 ENCODER (VIDEO)
INTERFACE	USB 2.0

DISTRIBUTED BY:

ECLERIS 