





## Guardian™ Electrodes for Long-Term Monitoring

Vermed's Guardian™ ECG monitoring electrodes are designed for use in ICU/CCU and other long-term applications. Guardian electrodes are available in soft cloth, breathable tape or occlusive foam. They also use low chloride liquid gel and hypoallergenic adhesives to reduce the possibility of irritation associated with long-term use. For quality readings in all of your critical monitoring situations choose Guardian from Vermed.

Vermed also carries a complete line of Lead Wires, Chart Paper and Defib Pads

- **Latex Free**
- **Gentle hypoallergenic adhesives**
- **Ag/AgCl sensor - low chloride, wet gel hydrates the skin and produces fast, clear traces**



Part Number	Description	Number of Electrodes			
		Per Strip	Per Pouch	Per Box	Per Case
A10006-3	2 1/4" Diameter Foam Metal Snap	3	3	75	600
A10006-60	2 1/4" Diameter Foam Metal Snap	3	60	-	600
A10006-3T	2 1/4" Diameter Tape Metal Snap	3	3	75	600
A10006-60T	2 1/4" Diameter Tape Metal Snap	3	60	-	600
A10007-3	2 1/4" Diameter Cloth Metal Snap	3	3	75	600
A10007-3NB	2 1/4" Diameter Cloth Metal Snap	3	3	-	600

To place an order or receive sample electrodes call your sales representative or Customer Service today!

Tenby Medical

A DIVISION OF PEKISKO HOLDINGS INC.

443 Village Square  
Sherwood Park, AB  
Canada T8A 4K2

Tel/Fax: 780-570-5755  
[tenby-medical.com](http://tenby-medical.com)