

EASY, TOUCH VITALS BY MIDMARK



IQ
vitals

MIDMARK®



Midmark IQvitals™ is an affordable, Windows® based, touchscreen device that provides automated blood pressure, temperature, SpO₂ and pulse rate readings with a touch of the screen. An alphanumeric keyboard allows for fast and easy entry of patient data, and the intuitive design increases productivity and simplifies staff training. IQvitals™ can be used as a stand-alone unit or will connect easily to a computer using IQmanager™ software.

TECHNOLOGICALLY ADVANCED, USER DESIGNED VITAL SIGNS

Efficiency

- Automated, all-in-one vital signs device saves time.
- Seamless connection to IQmanager™ software reduces transcription errors.
- Lightweight, portable device weighs just 3.9 lbs.
- Long battery life allows vitals to be captured continuously and quickly.
- Connects to select digital scales – use either digital or manual option to enter weight.
- Easily enter and view all vital signs data in one location, saving time.
- Manual vitals entry options for respiration rate, height, weight and pain score; calculates BMI automatically.

Technology

- Windows® based platform provides flexible, easy-to-upgrade software for the life of the product.
- Simplified touchscreen display for effortless, one-touch operation.
- Alphanumeric keyboard allows for quick and easy patient data entry.
- Intuitive user interface simplifies staff training.
- Meets or exceeds all standards for blood pressure, temperature, pulse and SpO₂.
- Built with hospital grade components to stand up to the demands of a busy practice.

Adaptability

- USB and serial connectivity options are standard.
- 100-reading memory for easy review, printing or download.
- Accessories available for patient testing, from pediatric to adult.
- Mobile cart, table top or wall mounting options available.
- External printer available.



Flexible mounting options to meet your needs.

IQvitals™ Technology by Midmark

Use as a standalone unit or connect to a computer using IQmanager™ software

Automated, all-in-one vital signs device saves time and reduces transcription errors.



Full-color, touchscreen display is quick and easy to operate.



Lightweight and portable at just 3.9 pounds with table-top, wall or mobile cart mounting options.



Connected flexibility with USB and serial ports standard.

Make your practice more efficient with an IQ Digital Diagnostic System

Workflow efficiency and improvements in patient care are just a few of the benefits that can result from implementing an IQ Digital Diagnostic System. Call us today and we'll put together a customized suite of solutions to help grow as your practice grows.

Diagnostic Products									
	ECG	Holter	Stress	Spirometer	Vital Signs	EHR/EMR Connectivity	Single Platform (IQmanager™)	Stand-alone	Thin - Client (IQpath™)
IQecg™	■					■	■		■
IQflex™	■							■	
IQclassic™	■							■	
IQholter™		■				■	■		
IQholter™ EX		■				■	■		
IQholter™ EP		■				■	■		
IQstress™			■			■	■		
IQspiro™				■		■	■		■
IQvitals™					■		■	■	
IQvitals™ PC					■		■		

System Requirements

- Windows® Vista Home, Windows® XP Pro, Windows® 7 Home Premium, 32-bit and 64-bit
- Intel® Core™ 2 Duo Processor E4300 or 64-bit processor or faster
- 2 GB of RAM
- 2 GB available hard disk space
- USB or Serial Port
- Mouse, keyboard, and high resolution display monitor

Please see the IQvitals™ Operation Manual or visit midmark.com for more detailed system requirements

Ordering Information

- 4-000-0500 IQvitals™ with Blood Pressure and Alaris® Turbo Temp® Thermometer
- 4-000-0510 IQvitals™ with Blood Pressure, Alaris® Turbo Temp® Thermometer and SpO₂
- 1-100-1605 IQvitals™ Thermal Printer
- 1-100-1603 Fairbanks® TeleWeigh™ Digital Scale
- 3-004-2000 IQvitals™ Mobile Cart

*Alaris and Turbo Temp are registered trademarks of Cardinal Health, Inc.
Windows is a registered trademark of Microsoft Corporation.
Intel Core is a trademark of Intel Corporation in the U.S. and/or other countries.
Fairbanks and TeleWeigh are trademarks of Fairbanks Scales, Inc.*

For information, pricing or to place an order, contact us at:
1-800-624-8950 or
1-310-516-6050 or
midmark.com

Midmark Diagnostics Group is an ISO 13485 Certified Company.

Midmark Diagnostics Group
1125 W. 190th Street, Gardena, CA 90248
midmark.com
1-800-624-8950

© 2010 Midmark Corporation
60 Vista Drive P.O. Box 286 Versailles, Ohio 45380-0286
Products subject to improvement changes without notice
Litho in U.S.A. 007-0836-00 Rev. A1 (8/10)



RETHINK THE CLINICAL SPACE

At the core of Midmark is a new way to think about caregiver/patient interaction and the delivery of efficient patient care. We are dedicated to working together with you to redefine the future of the clinical space.

midmarkclinicalsolutions.com

