


# Introducing Repositionable NeuroPlus

*Vermed's new line of repositionable disposable electrodes for use in Sleep Studies, EMG, NCS and EP.*

- NeuroPlus 43 is a small (2.5 cm square), aggressive electrode designed for easy placement and superior adhesion. Specially formulated adhesive solid gel ensures full sensor contact and excellent readings.
- NeuroPlus 40, 41 and 42 are designed with a more gentle adhesive gel for use in sensitive areas. Specially formulated adhesive gel is also repositionable to ensure proper placement. Sizes range from 2.2 cm square to 3.5 cm square.



A10040, A10041, and A10042, shown actual size

## Features

- Solid gel
- Small size
- Available in aggressive or gentle adhesive
- Hypoallergenic
- Cloth base
- Silver/Silver Chloride sensor
- Multipurpose

## Benefits

- Repositionable, easy clean up
- Easy placement on face or difficult locations
- Suitable for many different situations
- Long term use with no irritation
- Comfortable and conforming
- Clear signal
- Use for EMG or EKG



A10043, shown actual size



A10050-60, shown actual size

**All NeuroPlus electrodes meet or exceed AAMI specifications and are also suitable for use with EKG monitors.**

**NeuroPlus electrodes are latex free.**

Part Number	Description	Number of Electrodes		
		per strip	per pouch	per case
A10040-60	7/8" x 7/8" Cloth	3	60	600
A10041-60	1-3/8" x 7/8" Cloth	3	60	600
A10042-60	1-3/8" x 1-3/8" Cloth	3	60	600
A10043	1" x 1" Cloth	3	30	300
A10050-60	1" Diameter Foam Wet Gel	3	60	600

*To place an order or for more information, call your sales representative or customer service today!*

Tenby Medical

A DIVISION OF PEKISKO HOLDINGS INC.

443 Village Square  
Sherwood Park, AB  
Canada T8A 4K2

Tel/Fax: 780-570-5755  
[tenby-medical.com](http://tenby-medical.com)