

The Leader In
PFT Filters!


PulmoGUARD™

Fits most PFT devices without adapters

During a Pulmonary Function Test, bacteria and viruses can be expelled, posing serious risk to patients and health care workers alike. With a cross-contamination efficiency of 99.999%, PulmoGuard is safe, highly efficient...and very affordable.

PulmoGUARD™

INFECTION PROTECTION



- Meets ATS air resistance standards
- Dead space of 60cc
- Pleated media for better functionality
- Compact, non-intimidating design
- Transparent for easy, 100% visual inspection
- Fits most spirometers with no adapter

Please see back for Technical Specifications and uses.


PulmoGuard Part Numbers & Packaging

Part No.	Description	Qty	Packaging
29-7950-050	PulmoGuard Filter	50	Box
29-7950-100	PulmoGuard Filter	100	Box
29-7950-200	PulmoGuard Filter	200	4 Boxes/50
29-7950-300	PulmoGuard Filter	300	Box

Please see back for Kit options

PulmoGUARD II™

INFECTION PROTECTION



- All the benefits of the original PulmoGuard
PLUS –
- Fits many PFT instruments with no adapter
 - Dead space of 64cc

Please see back for Technical Specifications and uses.

PulmoGuard II Part Numbers & Packaging

Part No.	Description	Qty	Packaging
29-7920-050	PulmoGuard II Filter	50	Box
29-7920-100	PulmoGuard II Filter	100	Box
29-7920-200	PulmoGuard II Filter	200	4 Boxes/50
29-7920-300	PulmoGuard II Filter	300	Box

Please see back for Kit options

Your Dealer is:

Tenby Medical

A DIVISION OF PERSKO HOLDINGS INC.

1-888-309-0277

www.tenby-medical.com

sdi

DIAGNOSTICS

The Spirometry Source

Ten Hampden Drive, Easton MA 02375

www.sdiagnostics.com

PULMOGUARD TECHNICAL SPECIFICATIONS

PulmoGuard fits the following instruments without an adapter:
Adapters available for instruments not listed.

Assess - Peak Flow Meter
Burdick - Spirosense
Circadian - Models (if using black Tamarac pneumotach, need 7975)
Collins - Survey, Survey II, Survey III, Survey PLUS, Eagle, Eagle II
Cybermedics - Spinnaker TL, SL, CM3, CM310 (if white or black cone use 7980)
Cybermedics - Moose
Godart - All Models
Gould - 5004
Healthscan - Assess Peak Flow Meter
Jaeger - All Models except Spiro Pro
Jones - All Models except Satellite (small end of filter)
Koko - All Models
Med Graphics - 1070, 1085 (w/Prevent, use PG II; if no Prevent, 7970-spirometry, 7961-vol/diff)
Med Science - 570, 6257
Med Science - 3000, 4000 (spirometry-no adapt./vol-7961/diff-7960B)
Med Systems - Spirolite 101 & 201 (small end of filter), 343
New Med - 7000 I (Spirobyte)
PK Morgan - Benchmark (spirometry-7978, vol&diff-7960B), Flexi Flow (small end of filter), Trans Flow
Pneumedics - All Models (small end of filter)
Riko - AS300
S&M - CF-150B, GS520, Keystone III, Portascreen, Pulmoscreen
ViaSys/Sensormedics - Horizon, Ohio, 922, 2130, 2200, 2400, 2450, 2800, 6200, 2100 (spirometry-7961, vol-7960B, diff-7960B), Vmax (can use PGII w/7976 or 7977)
Spirolink - All Models
Spirometrics - LTE, 2500 LTE
Spiroscreen - All Models (small end of filter)
Spirotech - All Models
Tamarac - All Models except Presto Flash

PULMOGUARD II TECHNICAL SPECIFICATIONS

PulmoGuard II fits the following instruments without an adapter:
Adapters available for instruments not listed.

Breco - All Models
CDX - Spiro 110 (remove black fitting), Spiro 850
Cavitron/Donti - All Models
Clement-Clarke - Mini-Wright Peak Flow Meter
Collins - All DS series pulmonary labs with UBV (use 7982 adapter for hose), CPL, GS, SurveyTach, Eaglet, Hawk
Cosmed - All Models
Creative Biomedics - All Models
Futuremed - Discovery, Spirobank, Spirovision III
Medical Equip Designs - All Models (remove black fitting)
Med Graphics - preVent, PF/DX (needs to be reprogrammed if no filter before)
Micro Medical/Direct - All Models
MIR - All Models
Multispiro - All Models except PC (remove black fitting)
nSpiro - Wright Peak Flow Meter
Ohio - All Models
PK Morgan - 232
Puritan-Bennett - PS600, 900, 900L, VS400
Riko - All Models except AS300
S&M - Keystone
SDI - SBOoffice, SBG, spirotel, Spirolab, Spirolab II, Astra series with 7983 adapter
Smiths Medical - All models
Spirometrics - All Models except LTE, 2500 LTE
VacuMed - Discovery, Micro, MicroPlus, MicroGP/DL, MicroLoop, Pony, Spiroscreen, Wizard
ViaSys - MicroLoop, MicroLab
Vitalograph - All Models except Presto Flash
Wright Professional Peak Flow Meter

PulmoGuard Kits

Part No.	Description	Qty	Packaging
29-7951-050	PulmoGuard Filter + Cardboard Mouthpiece	50	Box
29-7951-100	PulmoGuard Filter + Cardboard Mouthpiece	100	Box
29-7953-040	PulmoGuard Filter + Cardboard Mouthpiece + The Klip	40	Box
29-7953-080	PulmoGuard Filter + Cardboard Mouthpiece + The Klip	80	Box
29-7954-050	PulmoGuard Filter + The Klip	50	Box
29-7954-100	PulmoGuard Filter + The Klip	100	Box
29-7957-040	PulmoGuard Filter + Comfit-D Disposable Mouthpiece	40	Box
29-7957-080	PulmoGuard Filter + Comfit-D Disposable Mouthpiece	80	Box
29-7958-040	PulmoGuard Filter + Comfit-D Disposable Mouthpiece + The Klip	40	Box
29-7958-080	PulmoGuard Filter + Comfit-D Disposable Mouthpiece + The Klip	80	Box

PulmoGuard II Kits

Part No.	Description	Qty	Packaging
29-7921-050	PulmoGuard II Filter + 1" O.D. Mouthpiece	50	Box
29-7921-100	PulmoGuard II Filter + 1" O.D. Mouthpiece	100	Box
29-7922-050	PulmoGuard II Filter + 1-1/16" O.D. Mouthpiece	50	Box
29-7922-100	PulmoGuard II Filter + 1-1/16" O.D. Mouthpiece	100	Box
29-7923-040	PulmoGuard II Filter + 1" O.D. Mouthpiece + The Klip	40	Box
29-7923-080	PulmoGuard II Filter + 1" O.D. Mouthpiece + The Klip	80	Box
29-7924-050	PulmoGuard II Filter + The Klip	50	Box
29-7924-100	PulmoGuard II Filter + The Klip	100	Box
29-7925-040	PulmoGuard II Filter + 1-1/16" O.D. Cardboard Mouthpiece + The Klip	40	Box
29-7925-080	PulmoGuard II Filter + 1-1/16" O.D. Cardboard Mouthpiece + The Klip	80	Box
29-7927-040	PulmoGuard II Filter + Comfit-D Disposable Mouthpiece	40	Box
29-7927-080	PulmoGuard II Filter + Comfit-D Disposable Mouthpiece	80	Box
29-7928-040	PulmoGuard Filter + Comfit-D Disposable Mouthpiece + The Klip	40	Box
29-7928-080	PulmoGuard Filter + Comfit-D Disposable Mouthpiece + The Klip	80	Box