



## PerformancePlus™ Electrodes for Stress and Holter Testing

Vermed is known throughout the industry as the leader in quality and performance for disposable diagnostic electrodes. Our PerformancePlus line continues that tradition. In critical diagnostic situations, use one of Vermed's PerformancePlus electrodes to insure quality readings, zero artifact and no repeat procedures. PerformancePlus electrodes are available in tape, cloth or foam in a variety of shapes and sizes to suit your needs. An assortment of no-waste packaging options are also available to fit your facility's usage.



A10005

- Ag/AgCl, high chloride, wet gel electrodes hydrate the skin and produce fast, clear tracings
- Aggressive adhesive keeps electrodes on the patient, no motion artifact or repeat testing
- Hypoallergenic adhesive
- Latex Free

Try our Stress Kit. 10 electrodes in one pouch, each designed for a specific location. Pouch includes placement diagram and application instructions.




A10006-5



A10006-5T



A10007-5



A10012-5S

Part Number	Description	Number of Electrodes			
		Per Strip	Per Pouch	Per Box	Per Case
A10005	1 1/2" diameter, Foam Metal Snap	5	5	100	600
A10005NB	1 1/2" diameter, Foam Metal Snap	5	5	-	600
A10005-1-60	1 1/2" diameter, Foam Metal Snap	1	60	-	600
A10005-7NB	1 1/2" diameter, Foam Metal Snap	7	7	-	700
A10005-60	1 1/2" diameter, Foam Metal Snap	5	60	-	600
A10006-5	1 7/8" diameter, Foam Metal Snap	5	5	100	600
A10006-1-60S	1 7/8" diameter, Foam Metal Snap	1	60	-	600
A10006-5T	2 1/4" diameter, Tape Metal Snap	5	5	100	600
A10006-1-60TS	2 1/4" diameter, Tape Metal Snap	1	60	-	600
A10007-5	2 1/4" diameter, Cloth Metal Snap	5	5	100	600
A10007-5NB	2 1/4" diameter, Cloth Metal Snap	5	5	-	600
A10012-5S	1 7/8" x 1 1/4", Foam Metal Snap	5	5	100	600
A10012-60S	1 7/8" x 1 1/4", Foam Metal Snap	5	60	-	600
A10015-10S	Stress Test Kit	50 Kits Per Case, 500 Electrodes			

To place an order or receive sample electrodes call your sales representative or Customer Service today!

**Tenby Medical**  
A DIVISION OF PEKISKO HOLDINGS INC.

443 Village Square  
Sherwood Park, AB  
Canada T8A 4K2  
Tel/Fax: 780-570-5755  
[tenby-medical.com](http://tenby-medical.com)