



## Regal™

### Diagnostic Resting ECG Tab Electrode

Vermed's Regal™ tab electrode for resting ECG provides superior adhesion and quality traces. Our silver/silver chloride coating and conductive adhesive gel provide an electrode that meets and exceeds the AAMI standard for disposable pre-gelled ECG electrodes.



A10023

Actual size 1 1/8" x 1 7/16"

- Proprietary adhesive solid gel for optimal adhesion
- Single use (eliminates cross contamination)
- Increased size for larger surface area contact
- Pliable substrate conforms to body
- Adhesive solid gel requires minimal clean-up
- Larger surface area contact reduces impedance and provides superior adhesion
- Repositionable
- Packaged in convenient re-sealable pouches



Part Number	Description	Number of Electrodes		
		Per Strip	Per Pouch	Per Case
A10023	Ag/AgCl Adhesive Solid Gel Tab	10	100	4000

To place an order or receive sample electrodes call your sales representative or Customer Service today!

Tenby Medical

A DIVISION OF PEKISKO HOLDINGS INC.

443 Village Square  
Sherwood Park, AB  
Canada T8A 4K2

Tel/Fax: 780-570-5755  
[tenby-medical.com](http://tenby-medical.com)