



# Gel Knit Removable Dressings

Riemann, Goodman, and Mannheim Gel Knit Nasal Dressings promote optimal surgical outcomes, control minor bleeding, and allow for easy removal with minimal discomfort to the patient

## Gel Knit Nasal Dressings



Riemann Nasal Dressing – 3.0cm/4.0cm



Goodman Nasal Dressing – 5.5cm



Mannheim Nasal Dressing – 8.0cm



Soak dressing to hydrate Gel Knit hydrocolloid fabric

## Gel Knit Postoperative Surgical Dressings

- Promotes optimal surgical outcomes
- Controls minor bleeding
- Prevents adhesions
- Allows for easy removal and patient comfort

Catalog No.	Epistaxis Solution	Size	Units Per Box
RR 300	Riemann	3.0cm	20
RR 400	Riemann	4.0cm	20
RR 500	Goodman	5.5cm	20
RR 800	Mannheim	8.0cm	20



**ArthroCare ENT**  
 680 Vaqueros Avenue  
 Sunnyvale, CA 94085-3523  
 main phone 408-736-0224  
 order entry **800-797-6520**  
 order entry fax 888-994-2782  
[www.arthrocareENT.com](http://www.arthrocareENT.com)

**ArthroCare Europe AB**  
 Skeppsbron 2  
 111 30 Stockholm  
 Sweden  
 phone +46 8 546 172 00  
 fax +46 8 546 172 39  
[info@arthrocare.se](mailto:info@arthrocare.se)

**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.

© 2010 ArthroCare Corporation. ArthroCare and Gel Knit are trademarks and/or registered trademarks of ArthroCare Corporation.