



# SilveRest™

**Adult Diagnostic Electrodes for Resting EKG and Echo sound**

- **Quality Construction**
- **Superior Performance**
- **Economical Testing**

Vermed SilveRest™ electrodes are the standard in diagnostic electrodes. Superior silver/silver chloride technology and quality construction produce an electrode that performs every time. Excellent adhesion and clear traces mean no repeat testing. SilveRest electrodes are available in various shapes, sizes and packaging configurations to suit your needs.



A10009



A10009F



A10013



A10049

- **Ag/AgCl, high chloride content, wet gel electrodes for fast, clear traces**
- **Self-adhering electrodes, no muscle artifact or repeat testing**
- **Available in economical, pliable, easy-to-peel-off vinyl tape or easy-placement, more adhesive, occlusive foam**
- **Latex Free**
- **Hypoallergenic adhesive**
- **Can be defibrillated**

| Part Number                 | Description       | Number of Electrodes |           |         |          |
|-----------------------------|-------------------|----------------------|-----------|---------|----------|
|                             |                   | Per Strip            | Per Pouch | Per Box | Per Case |
| <i>Wet Gel Electrodes</i>   |                   |                      |           |         |          |
| A10009                      | 1 3/8" Vinyl Tape | 10                   | 10        | 200     | 1200     |
| A10009-100                  | 1 3/8" Vinyl Tape | 10                   | 100       |         | 1000     |
| A10009-100F                 | 1 3/8" Foam       | 10                   | 100       | -       | 1000     |
| A10013-10SNB                | Rectangular Foam  | 5                    | 10        | -       | 1000     |
| A10013-100S                 | Rectangular Foam  | 5                    | 100       | -       | 1000     |
| <i>Solid Gel Electrodes</i> |                   |                      |           |         |          |
| A10049                      | Rectangular Foam  | 5                    | 50        | 500     | 2000     |

*To place an order or receive sample electrodes call your sales representative or Customer Service today!*

**Tenby Medical**

A DIVISION OF PEKISKO HOLDINGS INC.

443 Village Square  
Sherwood Park, AB  
Canada T8A 4K2

Tel/Fax: 780-570-5755  
**tenby-medical.com**