



Vermed SilverRest™

Diagnostic Resting EKG Tab Electrode

Vermed SilverRest™ tab electrodes for resting EKG provide excellent adhesion and quality traces. Our silver/silver chloride coating and conductive adhesive gel provide an electrode that meets or exceeds the AAMI standard for disposable pre-gelled EKG electrodes.

- **Conductive adhesive gel**
- **Specially designed strain relief tab**
- **Single use** (*eliminates cross contamination*)
- **Pliable substrate conforms to body**
- **Solid gel requires minimal clean-up**
- **Ultra-low impedance** (<100 OHMS) *promotes rapid baseline stabilization*
- **Repositionable**
- **Packaged in convenient, re-sealable pouches**



Specially designed strain relief tab is pliable, allowing electrode to remain secure.



Order Number	Description	Number of Electrodes		
		Per Strip	Per Pouch	Per Case
A10022	Ag/AgCl Adhesive Solid Gel Tab	10	100	4000
A10028	Ag/AgCl Adhesive Solid Gel Tab	10	100	4000

To place an order or for more information, call your sales representative or customer service today!