





**Tender-Trode™**  
*Electrodes for Pediatric and Neonatal Use*



## Tender Care for Tender Patients

Vermed's Tender-Trode™ line represents a complete selection of electrodes for infant and pediatric use. Each electrode is carefully designed with gentle, hypoallergenic adhesives and soft conforming materials to ensure comfort and effectiveness on young patients.

The TenderTrode line offers an electrode for every infant and pediatric application: NICU, Diagnostics, APNEA and PICU.



### Infant "Breathable" Electrodes

Our infant tape electrodes are designed for patient comfort. A special breathable tape with gentle adhesive and low chloride hydrogel provide a conforming and non-irritating electrode. Excellent for use in APNEA monitoring and all other pediatric applications.


### Neonatal Electrode

Our neonatal electrode is designed for use in NICU, Apnea monitoring and surgery. The carbon snap does not interfere with nuclear imaging systems or generate heat in MRI.


### Thermal Probe Covers

Our thermal probe covers use a delicate adhesive designed for tender infant skin. The highly reflective surface deflects heat from radiant warmers.






A10040-60



A10041-60



A10042-60



A10053-S

### Soft Cloth Electrodes

Adhesive solid gel electrodes with soft-cloth material are available in various sizes to suit your need. Suitable for a variety of pediatric applications.

### Pre-Wired Electrodes

We offer a complete line of pre-wired electrodes for use in NICU and other neonatal applications. Available in various sizes as well as metallic or radio-translucent lead wires.



Keeping Watch with Tender Care



Part Number	Description		Number of Electrodes		
			Per Strip	Per Pouch	Per Case
<b>Wet Gel Electrodes</b>					
A10000	1 5/8" Round	Tape	1	1	600
A10001	1 5/8" Round	Tape	3	3	600
A10000-60	1 5/8" Round	Tape	3	60	600
A10000-2-60	1 5/8" Round	Tape	2	60	600
A10046-60-3%	1 1/2" Round	Foam	3	60	600
A10047-60	1 7/8" Round	Tape	3	60	600
<b>Solid Gel Electrodes</b>					
A10045-60	1 3/4" x 7/8" Rect.	Foam	3	60	600
A10040-60	7/8" Square	Cloth	3	60	600
A10041-60	7/8" x 1 3/8" Rect.	Cloth	3	60	600
A10042-60	1 3/8" Square	Cloth	3	60	600
A10050-60	1" Round	Foam	3	60	600
<b>Tabs</b>					
A10048	1" Square		10	100	4000
<b>Pre-wired Electrodes</b>					
A10053-S	7/8" Square	Cloth .060	3	30	300
A10060-SRT	7/8" x 1 3/8" Rect.	Cloth .060	3	30	300
A10010-S	7/8" x 1 3/8" Rect.	Cloth .060	3	30	300
A10061-SRT	7/8" Square with RT lead	Cloth .060	3	30	300
<b>Probe Cover</b>					
A10008	Thermal probe cover	Foam	1	50	500

*To place an order or receive sample electrodes call your sales representative or Customer Service today!*



A DIVISION OF PEKISKO HOLDINGS INC.

443 Village Square  
 Sherwood Park, AB  
 Canada T8A 4K2

Tel/Fax: 780-570-5755

[tenby-medical.com](http://tenby-medical.com)